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



















20

Aesculap®

Triangular bracket bending pliers in accordance with KÖHLER

Symbols on product and packages

Symbol	Explanation
	Caution, general warning symbol Caution, see documentation supplied with the product

Intended use

The triangular bracket bending pliers in accordance with KÖHLER are used in orthodontics for bending spring-hard wire \varnothing 0.7 mm.

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by, or on order of a physician!

- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- ▶ Read, follow, and keep the instructions for use.
- ▶ Use the product only in accordance with its intended use, see Intended use.
- ▶ Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- ▶ Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- ▶ Replace any damaged components immediately with original spare parts.

Safe operation

Risk of injury and/or malfunction!
▶ Always carry out a function check prior to using the product.

Risk of injury when using the product beyond the field of view!
▶ Apply the product only under visual control.

Note
The inner side of the forceps contains guide grooves for intake of the wire. These milled slots allow exact bending of the triangle if the instructions in the diagram are followed and the wire is always gripped in the respective groove, see Fig. 1.

- ▶ Insert the wire, length 6 cm \varnothing 0.7 mm, spring-hard, into the diagonal groove in the forceps and close the forceps firmly, see Fig. 2.
- ▶ Guide the wire with your thumb firmly against the forceps.
- ▶ Turn the forceps against the wire and pressure of the thumb until a sharp kink occurs at an angle of 60°, see Fig. 3.
- ▶ Open the forceps.
- ▶ Lie the wire in the diagonal groove. When doing so, ensure that the first kinking remains in the horizontal plane and the forceps touches the opposite side, see Fig. 4.
- ▶ Close the forceps.
- ▶ Guide the wire with your thumb firmly against the forceps.
- ▶ Turn the forceps against the wire and the pressure of the thumb until a kink occurs at an angle of 60°, see Fig. 3. When doing so, guide the wire slightly upwards or downwards by approximately 2° from the horizontal plane, see Fig. 5.
- ▶ Open the forceps.
- ▶ Lie the wire in the diagonal groove.
- ▶ Remain in the horizontal plane with the triangular form that has developed. The opposite angle touches the forceps, see Fig. 6.
- ▶ When bending again, guide the wire quickly under or over the start of the wire. There is now an equilateral triangle. The start of the wire and the ongoing end are touching and are running parallel to one another, see Fig. 7 and see Fig. 8.

Bending the vertical part

- ▶ Insert the triangle into the forceps so that the wire engages exactly in the grooves, on the right or left hand side, see Fig. 9 and see Fig. 10.
- ▶ Close the forceps firmly.
- ▶ Quickly bend the end of the wire on the "triangle section" to a right angle, see Fig. 9 and see Fig. 10.
- ▶ Tilt the end of the wire by 10° to the tip of the triangular bracket that has already developed, see Fig. 11.

Aligning the triangle to the horizontal plane

- ▶ Insert half of the triangle into the designated milled slot of the forceps, see Fig. 12.
- ▶ Lie your thumb on the other half of the triangle.
- ▶ In order to parallelize, move thumb and forceps against each other, until the basal area of the triangle forms a plane. The vertical part of the wire is at a 90° angle to the basal area of the triangle, see Fig. 12 and see Fig. 13.
- ▶ Insert the triangle section into the triangle that has developed so that the vertical part in perpendicular milling and the buccal side of the later triangular bracket is in the diagonal groove, see Fig. 14.
- ▶ Close the forceps firmly.
- ▶ Bend the transversal part of the triangular section to the tip of the triangle until the basal area of the triangle and the transversal part run almost parallel.

Note
A slight divergence is good for later support of the triangular bracket on the teeth, see Fig. 15 and see Fig. 16.

- ▶ Ensure that the transversal part of the bracket runs exactly bisecting the angle over the tip of the triangle, see Fig. 17.

Note
The rigid part of a triangular bracket is, as a rule, distally against the physiological mesial drift of the teeth.

Note
Triangular bracket for the right and left side respectively of an upper or lower jaw, see Fig. 18 and see Fig. 19 respectively.
Finished triangular bracket, see Fig. 20.

Validated reprocessing procedure

General safety notes

Note
Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note
Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note
Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this. The recommended chemistry was used for validation.

Note
If there is no final sterilization, then a virucidal disinfectant must be used.

Note
For the latest information on reprocessing and material compatibility see also the Aesculap extranet at www.extranet.bbraun.com
The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.
Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- ▶ Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-1.org, link to Publications, Red Brochure – Proper maintenance of instruments.

Disassembling the product before carrying out the reprocessing procedure

- ▶ Disassemble the product immediately after use, as described in the respective instructions for use.
- ▶ Open up products with hinges.

Preparations at the place of use

- ▶ If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- ▶ Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- ▶ Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

Preparation before cleaning

- ▶ Disassemble the product prior to cleaning.

Cleaning/disinfection

Product-specific safety notes on the reprocessing procedure

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!
▶ Use cleaning and disinfecting agents according to the manufacturer's instructions which
– are approved for (e.g. aluminum, plastics, high-grade steel),
– do not attack softeners (e.g. in silicone).
▶ Observe specifications regarding concentration, temperature and exposure time.
▶ Do not exceed the maximum permissible cleaning temperature of 90 °C.

- ▶ Do not use oxidizing chemicals (e.g. H₂O₂), which could cause bleaching/layer loss of the product.
- ▶ Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.
- ▶ Mount jaws protection on the product.
- ▶ Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfecting.
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting.
 - as an integrated mechanical support measure for mechanical cleaning/disinfecting.
 - for additional cleaning of products with residues left after mechanical cleaning/disinfecting.
- ▶ Clean and disinfect microsurgical products mechanically if they can be placed securely in the machine or on the positioning aids.

Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Manual cleaning with immersion disinfection	<ul style="list-style-type: none"> Suitable cleaning brush 20 ml disposable syringe Keep working ends open for cleaning. When cleaning products with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning. Drying phase: Use a lint-free cloth or medical compressed air Drying phase: Use a lint-free cloth 	Chapter Manual cleaning/disinfection and subsection: <ul style="list-style-type: none"> Chapter Manual cleaning with immersion disinfection
Mechanical alkaline cleaning and thermal disinfection	<ul style="list-style-type: none"> Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). 	Chapter Mechanical cleaning/disinfection and subsection: <ul style="list-style-type: none"> Chapter Mechanical alkaline cleaning and thermal disinfection

Manual cleaning/disinfection

Risk to patients!
▶ The product must only be cleaned mechanically!

- Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- After manual cleaning/disinfection, check visible surfaces visually for residues.
- Repeat the cleaning /disinfection process if necessary.

Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water
FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
RT: Room temperature
*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Drain any remaining water fully.

Phase III

- Fully immerse the product in the disinfectant solution.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- Rinse/flush the product thoroughly (all accessible surfaces).
- Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- Rinse lumens with an appropriate disposable syringe at least five times.
- Drain any remaining water fully.

Phase V

- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

Mechanical cleaning/disinfecting

Note
The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note
The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	<ul style="list-style-type: none"> Concentrate, alkaline: <ul style="list-style-type: none"> pH = 13 <5 % anionic surfactant 0.5 % working solution <ul style="list-style-type: none"> pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water
FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection, maintenance and checks

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!
▶ Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- Dry the product if it is wet or damp.
- Repeat cleaning and disinfection of products that still show impurities or contamination.
- Check that the product functions correctly.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- This product should be assembled prior to sterilisation
- Check for compatibility with associated products.

Packaging

- Appropriately protect products with fine working tips.
- Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- Pack trays appropriately for the intended sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Steam sterilization

- Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- Validated sterilization process
 - Steam sterilization using fractional vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractional vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10⁻⁶, Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Prevacuum	270 °F/275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Technical Service

Risk of injury and/or malfunction!

► Do not modify the product.

► For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Modifications carried out on medical technical equipment may result in loss of warranty rights and applicable licenses.

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 (7461) 95-1602

Fax: +49 (7461) 16-5621

E-Mail: ats@aesculap.de

Or in the US:

Aesculap Inc.

Attn. Aesculap Technical Services

615 Lambert Pointe Drive

Hazelwood

MO, 63042

Aesculap Repair Hotline

Phone: +1 (800) 214-3392

Fax: +1 (314) 895-4420

Other service addresses can be obtained from the address indicated above.

Distributor in the US/Contact in Canada for product information and complaints

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