

Agility® Buccal Tubes Precision Manufactured By G&H® Orthodontics



Comfort is in the design.


The Agility® buccal tube is a miniaturized cast tube with rounded surfaces for exceptional patient comfort and acceptance. The design and manufacturing technologies ensure maximum strength for increased clinical efficiency and patient comfort. Both occlusal interference and gingival impingement are reduced, and hygiene is improved.

Single Non-Convertible

Maxillary 1st & 2nd Molar



Weldable



Bondable



Torque	Offset	Rx	M/D	R/L	Weldable		Bondable	
					.018	.022	.018	.022
0°	0°	-	3.5 mm	R	AS1-821-X	AS1-221-X	AS1B-821E-X	AS1B-221E-X
				L	AS1-822-X	AS1-222-X	AS1B-822E-X	AS1B-222E-X
-10°	8°	Roth ¹	3.5 mm	R	AS1-825-X	AS1-225-X	AS1B-825E-X	AS1B-225E-X
				L	AS1-826-X	AS1-232-X	AS1B-826E-X	AS1B-226E-X
-10°	0°	MBT ²	3.5 mm	R	AS1-831-X	AS1-231-X	-	-
				L	AS1-832-X	AS1-232-X	-	-
-14°	10°	MBT ²	3.5 mm	R	AS1-833-X	AS1-233-X	AS1B-833E-X	AS1B-233E-X
				L	AS1-834-X	AS1-234-X	AS1B-834E-X	AS1B-234E-X

Mandibular 1st & 2nd Molar

Torque	Offset	Rx	M/D	R/L	Weldable		Bondable	
					.018	.022	.018	.022
0°	0°	-	3.5 mm	R	AS1-822-X	AS1-222-X	AS1B-822E-X	AS1B-222E-X
				L	AS1-821-X	AS1-221-X	AS1B-821E-X	AS1B-221E-X
-10°	0°	-	3.5 mm	R	AS1-832-X	AS1-232-X	AS1B-832E-X	AS1B-232E-X
				L	AS1-831-X	AS1-231-X	AS1B-831E-X	AS1B-231E-X
-25°	5°	Roth ¹	3.5 mm	R	AS1-827-X	AS1-227-X	AS1B-827E-X	AS1B-227E-X
				L	AS1-828-X	AS1-228-X	AS1B-828E-X	AS1B-228E-X

Learn More -

Find out more about the Agility Tubes and the entire Agility family -

GHOrthodontics.com/Agility

¹G&H® prescriptions are not implied to be an exact version of any other system, nor do we claim any endorsement by the doctor.
²MBT is a trademark of 3M Unitek Corporation.