

## Seamlessly working together with our miniPrevail®TWIN and SL brackets

Our new miniPrevail<sup>®</sup>TUBE has been engineered with the patient and orthodontist in mind. Our distinct "crosshair" design takes the guesswork out of placement while the tube features offer easy wire placement and engagement. Finally, the low profile and sleek hook is designed for the utmost patient comfort.



GHOrthodontics.com TOLL-FREE 800-526-1026 or +1 317-346-6655 orders@ghortho.com

## **miniPrevail®TUBE**

Available in .018 & .022 - Single, Non-Convertible - Right & Left - Weldable & Bondable Roth\*,  $MBT^*$  or Standard prescriptions

1 <sup>st</sup> Molar Maxillary	Torque	Offset	D/M	₽	1 <sup>st</sup> Molar Mandibula	Torque	Offset	M/D	Q
Standard Rx	0°	0°	4.2 mm	•	Standard Rx	0°	0°	4.2 mm	•
Roth <sup>*</sup> Rx	-10°	8°	4.2 mm	•	Roth <sup>*</sup> Rx	-25°	5°	4.2 mm	•
MBT*	-14°	10°	4.2 mm	•	MBT*	-20°	0°	4.2 mm	•
MBT*	-10°	0°	4.2 mm	•	MBT*	-10°	0°	4.2 mm	•

2 <sup>nd</sup> Molar Maxillary	Torque	Offset	Q/W	₽
Standard Rx	0°	0°	3.2 mm	•
Roth* Rx	-10°	8°	3.2 mm	•
MBT*	-14°	10°	3.2 mm	•
MBT*	-10°	0°	3.2 mm	•

2 <sup>nd</sup> Molar Mandibular	Torque	Offset	Q/W	Q
Standard Rx	0°	0°	3.2 mm	٠
Roth* Rx	-25°	5°	3.2 mm	•
MBT*	-20°	0°	3.2 mm	•
MBT*	-10°	0°	3.2 mm	•

Tubes are available in packages of 10.

## Learn More -GHOrthodontics.com/miniPrevailTUBES

\*G&H<sup>®</sup> prescriptions are not implied to be an exact version of any other system, nor do we claim any endorsement by the doctor. MBT is a trademark of 3M Unitek Corporation.

Shop our full line of G&H manufactured products @ GHOrthodontics.com



BRACKETS | BANDS | TUBES | WIRES | SPRINGS | ELASTOMERICS